

Your Sportguard

requires minimal yet essential care

After removing your Sportguard, it's an excellent idea to clean it under warm, running tap water. Use a SOFT toothbrush and a non-abrasive toothpaste. Avoid using toothpaste indicated for whitening as the majority of them contain abrasive material that will damage your Sportguard. After cleaning, rinse the Sportguard thoroughly, allow it to air dry and store it in the container provided without soaking.

Important Notes:

- Store your Sportguard out of reach of your family pets; they find these appliances very difficult to resist chewing once their owner has used it.
- Do not chew on your Sportguard because you could weaken the effectiveness of the material.
- Do not share your Sportguard with anyone.
- Avoid high temperatures and direct sunlight to minimize distortion.



OTTAWA 67'S
OFFICIAL SPORTSGUARD MANUFACTURER

- If the Sportguard absorbs a sharp blow, it may need replacement. It may no longer be as protective as it was, even if it still looks intact.
- The Sportguard material deteriorates over time and loses resilience, reducing its effectiveness. Therefore, they require replacement every two-three years. Replace earlier if it becomes cracked, torn, split, or significant wear appears, or there is poor retention.
- If you have any additional questions or concerns, be sure to bring them to the attention of your dental care providers.